

L-Arginyl-L-arginine 7-Amido-4-methylcoumarin Trihydrochloride

Art. ID TRC-A769523-250MG

Unit 250 mg

Description

Application Notes: L-Arginyl-L-arginine 7-Amido-4-methylcoumarin Trihydrochloride is a fluorogenic substrate for cathepsin B or aminopeptidase III. It yields a blue fluorescent solution upon cleavage.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Arginyl-L-arginine 7- Amido-4-methylcoumarin Trihydrochloride	[201847-69-0]				