

8-Azidoadenosine 3':5'-Cyclic Monophosphate

Art. ID TRC-A822003-25MG
Unit 25 mg

Description

Category: Photoaffinity Labels /// Application Notes: 8-Azidoadenosine 3':5'-Cyclic Monophosphate is a photoaffinity label for cAMP-dependent enzymes. Specifically it has shown its utility as a photoaffinity probe for cAMP binding sites in human erythrocyte membranes.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	8-Azidoadenosine 3':5'- Cyclic Monophosphate	[31966-52-6]				