

N-Acetyl-3,3',5-triiodo-L-thyronine-d3

Art. ID TRC-A188887-25MG
Unit 25 mg

Description

Category: Isotope Labelled Compounds, Aromatics, Amino Acids and Derivatives, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: N-Acetyl-3,3',5-triiodo-L-thyronine-d3 is the isotope labelled analog of N-Acetyl-3,3',5-triiodo-L-thyronine. N-Acetyl-3,3',5-triiodo-L-thyronine is used in biological studies to evaluate the thyroxine deiodinase of liver microsomes inhibition in rat models. /// References: Koehrie, J., et al.: J. Biol. Chem., 261, 11613 (1986)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-3,3',5-triiodo -L-thyronine-d3					