

L-Alanyl-L-proline-d3 t-Butyl Ester Hydrochloride

Art. ID TRC-A511347-25MG
Unit 25 mg

Description

Application Note: L-Alanyl-L-proline-d3 t-Butyl Ester Hydrochloride is an intermediate in synthesizing Enalapril-d3 Hydrochloride (E555258), which is the hydrochloride salt of Enalapril-d3 (E555253), a labeled version of Enalapril (E555403), an antihypertensive that is also an angiotensin-converting enzyme (ACE) inhibitor.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Alanyl-L-proline-d3 t -Butyl Ester Hydrochlor ide					