

(4-Aminophenyl)(piperazin-1-yl)methanone

Art. ID TRC-A625900-750MG

Unit 750 mg

Description

Application Notes: (4-Aminophenyl)(piperazin-1-yl)methanone is used as a reagent in the synthesis of novel (triazolyl)(2,4-difluorophenyl)(4-substituted-1-piperazinyl)-2-propanol derivatives which can exhibit antifungal activity. /// References: Sheng, C., et al.: Yaoxue Xuebao, 39, 984 (2004)| Sheng, C., et al.: Chinese Chem. Lett., 15, 404 (2004)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(4-Aminophenyl)(piperaz in-1-yl)methanone	[72141-42-5]				