

(3S,4S)-4-Aminotetrahydro-3-furanol

Art. ID TRC-A629715-500MG
Unit 500 mg

Description

Category: Amines, Heterocycles, Intermediates /// Application Notes: (3S,4S)-4-Aminotetrahydro-3-furanol, is a novel, optically active heterocyclic amino alcohols, that can be used as building blocks for biologically active molecules. These compounds, and their derivatives can also serve as both monosaccharide mimetic structures, and amino acid surrogates. /// References: Moon Kim, B., et al.: tetra. Lett., 39, 6921 (1998)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(3S,4S)-4-Aminotetrahydro-3-furanol	[214629-29-5]				