

5-Azido Uridine

Art. ID TRC-A926800-5MG
Unit 5 mg

Description

Category: Nucleotides, Bases and Related Reagents, Pharmaceuticals, Intermediates and Fine Chemicals ///
Appearance: Yellow Solid /// Application Notes: Used in a method and kit for labeling nucleic acid in living cells. The labeled nucleic acid can be used for detecting cell proliferation and cell apoptosis, labeling plasmid vector, virus, small interfering RNA and microRNA, and studying early development.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	5-Azido Uridine	[1355028-82-8]				