

Adenosine 3',5'-Bisphosphate Dicalcium Hydrate

Art. ID TRC-A280490-250MG

Unit 250 mg

Description

Category: Amines, Nucleotides, Bases and Related Reagents, Phosphorylating and Phosphitylating Agents ///

Application Notes: Adenosine 3',5'-Bisphosphate Dicalcium Hydrate is a substrate for the amidohydrolase superfamily. ///
References: Cummings, J.A., et. al.: Biochemistry, 53, 591 (2014)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adenosine 3',5'-Bisphosphate Dicalcium Hydrate					