

### **(R)-Ambrisentan-d3**

Art. ID TRC-A575874-100MG  
Unit 100 mg

#### Description

Application Notes: (R)-Ambrisentan-d3 is the deuterated version of (R)-Ambrisentan (A575873), which is the R-isomer impurity of Ambrisentan (A575860) that is a nonpeptide endothelin ETA receptor antagonist, used as antihypertensive drug. // References: Riechers, H., et al.: J. Med. Chem., 39, 2123 (1996)| Billman, G.E., et al.: Curr. Opin. Invest. Drugs, 3, 1483 (2002)| Vatter, H., et al.: Clin. Neuropharmacol., 26, 73 (2003)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(R)-Ambrisentan-d3					