

4-Amino-2,6-dimethyl-8-(2'-deoxy-beta-D-ribofuranosyl)-7(8H)-pteridone

Art. ID TRC-A605170-100MG
Unit 100 mg

Description

Category: Fluorescent Labels and Indicators, Nucleotides, Bases and Related Reagents, /// Appearance: Crystalline Solid /// Application Notes: A highly fluorescent adenosine analogue, which in a dimethoxytrityl, phosphoramidite protected form, can be site-specifically inserted into oligonucleotides through a 3',5'-phosphodiester linkage using an automated DNA synthesizer. /// References: Hawkins, M.E., et al.: Analytical Biochem., 298, 231 (2001)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Amino-2,6-dimethyl-8-(2'-deoxy-beta-D-ribofuranosyl)-7(8H)-pteridone	[195442-56-9]				