

[[[3-[(2-Amino-1,6-dihydro-6-oxo-9H-purin-9-yl)methyl]phenyl]methoxy]methyl]-phosphonic Acid Monoammonium Salt

Art. ID TRC-A610075-50MG

Unit 50 mg

Description

Application Note: [[[3-[(2-Amino-1,6-dihydro-6-oxo-9H-purin-9-yl)methyl]phenyl]methoxy]methyl]-phosphonic Acid Monoammonium Salt is derived from Diisopropyl Phosphite (D459370), which is used as antiwear lubricant additive together with triphenyl thiophosphate on T8 steel/Al₂O₃ ceramics.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	[[[3-[(2-Amino-1,6-dihydro-6-oxo-9H-purin-9-yl)methyl]phenyl]methoxy]methyl]-phosphonic Acid Monoammonium Salt	[143320-22-3]				