

## **ASPT**

Art. ID TRC-A780108-2.5G  
Unit 2.5 g

### Description

Application Notes: ASPT is an (aminostyryl)pyridinium hemicyanine dye that has solvent-dependent one- and two-photon emissions. ASPT has been utilized in the synthesis of fluorescent poly(?-aprolactone), where the alcohol group serves as the initiator of the reaction.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ASPT	[159721-38-7]				