

(2S,3R)-2-amino-3-hydroxy-4-(4-nitrophenyl)butanoate

Art. ID TRC-A938100-250MG

Unit 250 mg

Deliverydetails to be checked

Description

Application Note: (2S,3R)-2-amino-3-hydroxy-4-(4-nitrophenyl)butanoate is a useful chiral building block, primarily used in the synthesis of beta-lactone antibiotics.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2S,3R)-2-amino-3-hydroxy-4-(4-nitrophenyl)butanoate	[99268-38-9]				