

3-[(6-Acetyl-2-naphthalenyl)amino]alanine

Art. ID TRC-A186960-500MG
Unit 500 mg

Description

Category: Amines, Aromatics, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes:
3-[(6-Acetyl-2-naphthalenyl)amino]alanine is used in the synthesis of genetically encodable fluorescent unnatural amino acid L-3-(6-Acetylnaphthalen-2-ylamino)-2-aminopropanoic acid. /// References: Xiang, Z., et al.: J. Org. Chem., 76, 6367 (2011)| Lee, H.S., et al.: J. Am. Chem. Soc., 131, 12921 (2009)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-[(6-Acetyl-2-naphthalenyl)amino]alanine	[1185251-08-4]				