

3-Acetylpyridine-5-boronic Acid Pinacol Ester

Art. ID TRC-A188305-2.5G
Unit 2.5 g

Description

Application Notes: 3-Acetylpyridine-5-boronic Acid Pinacol Ester is a boron containing reagent used in the synthesis of pyrimidine derivatives as inhibitors of 6-phosphofructose-2-kinase used in the treatment of cancer.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Acetylpyridine-5-boronic Acid Pinacol Ester	[1103862-13-0]				