

N-Acetyl-3,3',5'-triiodo-L-thyronine-d3

Art. ID TRC-A188882-25MG
Unit 25 mg

Description

Category: Isotope Labelled Compounds, Aromatics, Amino Acids and Derivatives, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: N-Acetyl-3,3',5'-triiodo-L-thyronine-d3 is the isotope labelled analog of N-Acetyl-3,3',5'-triiodo-L-thyronine. N-Acetyl-3,3',5'-triiodo-L-thyronine is used in biological studies for iodothyronine deiodinase of liver inhibition /// References: Mol, J.A., et al.: Biochemical and Biophysical Res. Comm., 124, 475 (1984)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-3,3',5'-triiodo-L-thyronine-d3					