

DL-Alanine Methylamide-d3

Art. ID TRC-A510592-25MG
Unit 25 mg

Description

Category: Isotope Labelled Compounds, Amino Acids and Derivatives /// Application Notes: DL-Alanine Methylamide-d3 is the isotope analog of DL-Alanine Methylamide. DL-Alanine Methylamide is a reagent used in the acetylation of DL-amino acid derivatives, which are potent new agents for the treatment of epilepsy. /// References: Conley, J.D., et al.: J. Med. Chem., 30, 567 (1987)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-Alanine Methylamide-d3					