

(S)-2-Aminobut-3-en-1-ol Hydrochloride

Art. ID TRC-A602130-250MG

Unit 250 mg

Description

Application Notes: (S)-2-Aminobut-3-en-1-ol Hydrochloride is used as a reagent in the synthesis of bicyclic bisarylimidazole JNK3 inhibitors. It is also used in the preparation of kalkitoxin and its stereoisomers which are potent neurotoxins. /// References: Rech, J.C., et al.: J. Am. Chem. Soc., 129, 490 (2007)| Yokokawa, F., et al.: Tetrahedron, 60, 6859 (2004)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(S)-2-Aminobut-3-en-1-ol Hydrochloride	[219803-57-3]				