

1,3-Bis(allylamino)-2-propanol Dihydrochloride

Art. ID TRC-A614293-500MG
Unit 500 mg

Description

Application Note: 1,3-Bis(allylamino)-2-propanol Dihydrochloride is used in preparation of polymers, such as ethylenic polymers and crosslinked poly(allylamine) polymers.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,3-Bis(allylamino)-2-p ropanol Dihydrochloride	[1101184-94-4]				