

N-(3-Aminopropanoyl) Lenalidomide Hydrochloride

Art. ID TRC-A628105-100MG

Unit 100 mg

Description

Application Notes: N-(3-Aminopropanoyl) Lenalidomide is a derivative of Lenalidomide (L328000), an immunomodulatory drug, analog of Thalidomide. /// References: Bartlett, J.B., et al.: Nature Rev., 4, 314 (2004), Mitsiades, C.S., et al.: Curr. Opin. Invest. Drugs, 5, 635 (2004), List, A., et al.: N. Engl. J. Med., 352, 549 (2005)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-(3-Aminopropanoyl) Lenalidomide Hydrochloride					
	e					