

(2R)-2-Amino-3-(pyrrolidin-3-yl)propanoic Acid

Art. ID TRC-A629305-1G

Unit 1 g

Description

Application Notes: (2R)-2-Amino-3-(pyrrolidin-3-yl)propanoic Acid is a synthetic block made from (R)-2-Pyrrolidinone-5-carboxylic Acid (P998500), which is versatile reactant used in the synthesis of diphthamide. ///
References: Evans, D.A. and Lundy, K.M.: J. Am. Chem. Soc., 114, 1495 (1992)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2R)-2-Amino-3-(pyrroli din-3-yl)propanoic Acid	[1313183-08-2]				