

L-Alanine 7-Amido-4-methylcoumarin, Trifluoroacetate Salt

Art. ID TRC-A637323-2.5G
Unit 2.5 g
Deliverydetails /not restricted

Description

Application Notes: L-Alanine 7-Amido-4-methylcoumarin, Trifluoroacetate Salt is a fluorogenic substrate for aminopeptidase which yields a blue fluorescent solution. /// References: D. Mantle, et al.: Biochem. J. 1983, 211, 567

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Alanine 7-amido-4-methylcoumarin trifluoroacetate salt	[96594-10-4]				