

Androstane-3alpha,17beta-diol 17-Glucuronide-d5

Art. ID TRC-A637907-25MG
Unit 25 mg

Description

Application Notes: Androstane-3alpha,17beta-diol 17-Glucuronide-d5 is the isotope labelled analog of Androstane-3alpha,17beta-diol 17-Glucuronide. Androstane-3alpha,17beta-diol 17-Glucuronide is a derivative of androsterone (A637535), which is a steroid hormone with weak androgenic activity. It is metabolized from Testosterone (T155000) in the liver and is used in doping analysis to detect testosterone misuse. /// References: Park, J., et al.: J. Anal. Toxicol., 14, 66 (1990)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Androstane-3alpha,17beta-diol 17-Glucuronide-d5	[384380-31-8]				