

Anthracene 9,10-dipropionate Disodium Salt

Art. ID TRC-A678460-250MG

Unit 250 mg

Description

Application Notes: Anthracene 9,10-dipropionate Disodium Salt is Anthracene (A678400) based derivative. Anthracene 9,10-dipropionate Disodium Salt is used as reagent for determination of quantum efficiency of singlet oxygen generation in aqueous media. /// References: Kuznetsova, A., et al.: Russ. J. Gen. Chem., 71, 36 (2001)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Anthracene 9,10-dipropionate Disodium Salt	[82767-90-6]				