

AsCy3_E

Art. ID TRC-A780195-50MG
Unit 50 mg
Deliverydetails to be checked

Description

Application Note: AsCy3_E, is a cell permeable, biasenical fluorescent reagent that enables chemistry-driven approach for site-specific labeling of recombinant proteins in living cells avoiding the need for larger protein tags that can disrupt the structure or function of target proteins.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	AsCy3_E	[1421826-76-7]				