

4,4'-Azodianiline

Art. ID TRC-A940000-250MG
Unit 250 mg

Description

Category: Aromatics, Miscellaneous Reagents /// Appearance: Brown Solid /// Application Notes: 4,4'-Azodianiline is used in cross-linking studies of epoxy resins and the synthesis of surface relief holographic materials. ///

References: Goldenberg, L.M. et al.: J. Mat. Chem., 20, 9161 (2010) Yu, W.C. et al.: Macromol., 21, 365 (1988)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,4'-Azodianiline	[538-41-0]				