

Custom Solution - 10 components in HNO3 5%

Art. ID LM24-CUS-72964

Unit 1 x 100 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	10	mg/L		
	Iron (Fe)	[7439-89-6]	200	mg/L		
	Manganese (Mn)	[7439-96-5]	50	mg/L		
	Uranium (U)	[7440-61-1]	30	mg/L		
	Cadmium (Cd)	[7440-43-9]	5	mg/L		
	Vanadium (V)	[7440-62-2]	140	mg/L		
	Chromium (Cr)	[7440-47-3]	25	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Nickel (Ni)	[7440-02-0]	20	mg/L		
	Aluminium (Al)	[7429-90-5]	200	mg/L		