

Merck ICP multi-element standard V - 26 components in HCl 5%

Art. ID LM24-OP-4000034-500ML
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	20	mg/L		
	Arsenic (As)	[7440-38-2]	20	mg/L		
	Boron (B)	[7440-42-8]	2	mg/L		
	Barium (Ba)	[7440-39-3]	2	mg/L		
	Beryllium (Be)	[7440-41-7]	1	mg/L		
	Calcium (Ca)	[7440-70-2]	10	mg/L		
	Cadmium (Cd)	[7440-43-9]	2	mg/L		
	Chromium (Cr)	[7440-47-3]	2	mg/L		
	Copper (Cu)	[7440-50-8]	2	mg/L		
	Iron (Fe)	[7439-89-6]	2	mg/L		
	Mercury (Hg)	[7439-97-6]	5	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Lithium (Li)	[7439-93-2]	2	mg/L		
	Magnesium (Mg)	[7439-95-4]	1	mg/L		
	Manganese (Mn)	[7439-96-5]	1	mg/L		
	Sodium (Na)	[7440-23-5]	20	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Phosphorus (P)	[7723-14-0]	10	mg/L		
	Lead (Pb)	[7439-92-1]	20	mg/L		
	Scandium (Sc)	[7440-20-2]	1	mg/L		
	Selenium (Se)	[7782-49-2]	20	mg/L		
	Strontium (Sr)	[7440-24-6]	1	mg/L		
	Tellurium (Te)	[13494-80-9]	20	mg/L		
	Titanium (Ti)	[7440-32-6]	2	mg/L		
	Yttrium (Y)	[7440-65-5]	1	mg/L		
	Zinc (Zn)	[7440-66-6]	2	mg/L		