

Custom Solution - 13 components in HNO3 5%

Art. ID LM24-CUS-72626-2X250ML
Unit 2 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	2500	mg/L		
	Beryllium (Be)	[7440-41-7]	20	mg/L		
	Calcium (Ca)	[7440-70-2]	5000	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	2000	mg/L		
	Manganese (Mn)	[7439-96-5]	200	mg/L		
	Sodium (Na)	[7440-23-5]	4000	mg/L		
	Nickel (Ni)	[7440-02-0]	250	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	250	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		