

### **Silver (Ag) 1000 mg/kg in Base oil 20 cSt**

Art. ID                      LM24-PE-2000516  
Unit                         50 mL  
Deliverydetails            No Dangerous Good /not restricted

#### **Description**

Validity: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	1000	mg/kg		
	Silver (Ag)	[7440-22-4]	0,1	%		