

Beryllium (Be) 5000 mg/kg in Base oil 20 cSt

Art. ID LM24-PE-2000526
Unit 50 mL
Deliverydetails No Dangerous Good /not restricted

Description

Validity: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	5000	mg/kg		
	Beryllium (Be)	[7440-41-7]	0,5	%		