

## **Nickel (Ni) 5000 mg/kg in Base oil 20 cSt**

Art. ID                      LM24-PE-2000561  
Unit                         50 mL  
Deliverydetails            No Dangerous Good /not restricted

### Description

Validity: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nickel (Ni)	[7440-02-0]	5000	mg/kg		
	Nickel (Ni)	[7440-02-0]	0,5	%		