

Silicon (Si) 1000 mg/kg in Base oil 20 cSt

Art. ID LM24-PE-2000569
Unit 50 mL
Deliverydetails No Dangerous Good /not restricted

Description

Validity: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	1000	mg/kg		
	Silicon (Si)	[7440-21-3]	0,1	%		