

## **Tungsten (W) 1000 mg/kg in Base oil 20 cSt**

Art. ID                      LM24-PE-2000581  
Unit                         50 mL  
Deliverydetails            No Dangerous Good /not restricted

### Description

Validity: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tungsten (W)	[7440-33-7]	1000	mg/kg		
	Tungsten (W)	[7440-33-7]	0,1	%		