

Tungsten (W) 5000 mg/kg in Base oil 20 cSt

Art. ID LM24-PE-2000582
Unit 50 mL
Deliverydetails No Dangerous Good /not restricted

Description

Validity: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tungsten (W)	[7440-33-7]	5000	mg/kg		
	Tungsten (W)	[7440-33-7]	0,5	%		