

Custom Solution - 3 components 1000 mg/L each of in HNO3 5% (Part 2 out of 2)

Art. ID LM24-CUS-72955-2
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	1000	mg/L		
	Zirconium (Zr)	[7440-67-7]	1000	mg/L		
	Niobium (Nb)	[7440-03-1]	1000	mg/L		