

## **Tinned sausage (lyoner sausage)**

Art. ID MUVA-BR-1009  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

Best before:

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fat		22,64 ± 0,9	g/100g		
	Water (103 °C)		61,57 ± 0,59	g/100g		
	Protein N x 6,38 (Kjeldahl)		12,73 ± 0,25	g/100g		
	Chloride (Cl)		1121 ± 46	mg/100g		
Sodium chloride (Cl) calculated	Chloride (Cl)		1,85 ± 0,08	g/100g		
550 °C - 600 °C	Ash		2,55 ± 0,08	g/100g		
	Nitrate		21,38 ± 0,77	mg/kg		
Sodium (Na) calculated	Sodium (Na)	[7440-23-5]	1,90 ± 0,12	g/100g		
	Sodium (Na)	[7440-23-5]	746 ± 47	mg/100g	(AAS/ICP)	
	Hydroxyproline	[51-35-4]	0,2 ± 0,04	g/100g	photom.	