

Custom Solution -Merck ICP Multi-element standard Solution (XII) - 8 components 1000 mg/L each of in HCl 5%

Art. ID LM24-CUS-72975-4X30ML
Unit 4 x 30 mL bottle
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	1000	mg/L		
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Sulfur (S)	[7704-34-9]	1000	mg/L		
	Silicon (Si)	[7440-21-3]	1000	mg/L		
	Tungsten (W)	[7440-33-7]	1000	mg/L		
	Vanadium (V)	[7440-62-2]	1000	mg/L		
	Zirconium (Zr)	[7440-67-7]	1000	mg/L		