

Low alloy Steel for solid spectrometry, combustion and wet-way methods, Chips

Art. ID CZ-SPL-300A-Chips
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	Base			
	Carbon (C)	[7440-44-0]	0,0033	%		
	Manganese (Mn)	[7439-96-5]	0,039	%		
	Silicon (Si)	[7440-21-3]	0,0032	%		
	Phosphorus (P)	[7723-14-0]	0,0035	%		
	Sulfur (S)	[7704-34-9]	0,0021	%		
	Copper (Cu)	[7440-50-8]	0,008	%		
	Chromium (Cr)	[7440-47-3]	0,0092	%		
	Nickel (Ni)	[7440-02-0]	0,0075	%		
	Aluminium (Al)	[7429-90-5]	0,0017	%		
	Molybdenum (Mo)	[7439-98-7]	0,0015	%		
	Tungsten (W)	[7440-33-7]	~0,0017	%		
	Vanadium (V)	[7440-62-2]	~0,0011	%		
	Titanium (Ti)	[7440-32-6]	~0,0004	%		
	Cobalt (Co)	[7440-48-4]	0,0022	%		
	Arsenic (As)	[7440-38-2]	0,0027	%		
	Tin (Sn)	[7440-31-5]	~0,0008	%		
	Boron (B)	[7440-42-8]	~0,0002	%		
	Antimony (Sb)	[7440-36-0]	~0,0015	%		
	Lead (Pb)	[7439-92-1]	~0,0007	%		
	Zinc (Zn)	[7440-66-6]	~0,0007	%		
	Nitrogen (N)	[7727-37-9]		%		
	Bismuth (Bi)	[7440-69-9]	0,0028	%		
	Tantalum (Ta)	[7440-25-7]	~0,003	%		