

### **Custom Solution - 4 Components 100 mg/L each of in Methanol**

Art. ID LM24-CUS-72986  
Unit 5 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Atrazine D5	[163165-75-1]	100	mg/L		
	2,6-Dichlorobenzamide-3 ,4,5-d3	[1219804-28-0]	100	mg/L		
	Molinate D5	[n/a]	100	mg/L		
	Alachlor D13 (2,6-dieth ylphenyl D13)	[1015856-63-9]	100	mg/L		