

Custom Solution - 4 components in HNO3 1%

Art. ID LM24-CUS-72988
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	2000	mg/L		
	Chloride (Cl-)		3000	mg/L		
	SO42-		1000	mg/L		
	Phosphate (PO43-)		250	mg/L		