

Custom Solution - ICP Calibration Standard - 10 components in 2% HNO₃

Art. ID LM24-CUS-70135-3X500ML
Unit 3 x 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	10	mg/L		
	Potassium (K)	[7440-09-7]	5	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Aluminium (Al)	[7429-90-5]	1	mg/L		
	Copper (Cu)	[7440-50-8]	1	mg/L		
	Manganese (Mn)	[7439-96-5]	1	mg/L		
	Barium (Ba)	[7440-39-3]	0,2	mg/L		
	Calcium (Ca)	[7440-70-2]	0,2	mg/L		
	Magnesium (Mg)	[7439-95-4]	0,2	mg/L		
	Zinc (Zn)	[7440-66-6]	0,2	mg/L		