

Custom Solution - 5 components 1000 mg/L in HNO3 5%

Art. ID LM24-CUS-73001

Unit 100 mL bottle

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Tungsten (W)	[7440-33-7]	1000	mg/L		
	Zirconium (Zr)	[7440-67-7]	1000	mg/L		