

Custom Solution - 3 components 50 µg/mL each of in HNO3 2%

Art. ID LM24-CUS-70973-2X50ML
Unit 2 x 50 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	50	µg/mL		
	Gallium (Ga)	[7440-55-3]	50	µg/mL		
	Zirconium (Zr)	[7440-67-7]	50	µg/mL		