

**Custom Solution - Si 0,025 µg/mL in HNO3 2%**

Art. ID LM24-CUS-73006  
Unit 10 x 30 mL bottle  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

**Description**

Conservation 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,025	µg/mL		