

Hexamethylacetone 1000 µg/mL in Methanol

Art. ID LM24-N-13600-0999-1000ME5

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2,4,4-Tetramethyl-3-pentanone	[815-24-7]	1000	µg/mL		