

Tridecyl acrylate (mixed isomers) 1000 µg/mL in Methanol

Art. ID LM24-N-13710-0215-1000ME5

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tridecyl acrylate (mixed isomers)	[3076-04-8]	1000	µg/mL		