

EasyMP™ fragments (10-100 µm) in ethanol (95 vol%.) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-001
Unit 30 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1170 III E1

Description

available colors: Clear/Red/Blue/Violet/Fluorescent green (Violet/Green only for PE)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PP (polypropylene) / PVC (polyvinyl chloride) / PET (polyethylene terephthalate) / PS (polystyrene) / PUR (polyurethane) / TPU (thermoplastic polyurethane) / PS (polystyrene) / PA6 (Polyamide 6) / PA12 (Polyamide 12) / PC (polycarbonate) / CA (cellulose acetate) / PLA (polylactic acid) / PTFE (polytetrafluoroethylene) / POM (polyoxymethylene)		<1	%		