

Dry, neat (powder) fragments (500-1000 µm) - different polymer types and colors available - please contact us for details

Art. ID MPS-009
Unit 0,5 g

Description

available colors: Clear/Black/Red/Blue

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PP (polypropylene) / PVC (polyvinyl chloride) / P ET (polyethylene terephthalate) / PS (polystyrene) / TPU (thermoplastic polyurethane) / ABS (acrylonitrile butadiene styrene) / PA6 (Polyamide 6) / PA12 (Polyamide 12) / PMMA (Polymethyl methacrylate) / PC (polycarbonate) / CA (cellulose acetate) / PLA (polylactic acid) / PTFE (polytetrafluoroethylene) / POM (polyoxymethylene)					