

EasyMP™ beads (10-30 µm) in ethanol (95 vol%.) solution - different polymer types and colors available - please contact us for details

Art. ID MPS-012
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1170 III E1

Description

available colors: Clear/Red/Fluorescent green (Color only for PE)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
available polymer types	PE (polyethylene) / PS (polystyrene) / PMMA (P olymethyl methacrylate)		<1	%		